

# Patient Information Leaflet

## RESP-EASE<sup>TM</sup>

Sterile **6%** hypertonic  
sodium chloride solution.



**RESP-EASE<sup>TM</sup>**



**VENTURE**  
HEALTHCARE<sup>LTD</sup>



0459

Distributed by **Venture Healthcare Ltd**

Office: 16 Southborough Road, London, KT6 6JN

0345 468 0072 | [contact@venhealth.co.uk](mailto:contact@venhealth.co.uk) | [www.venhealth.co.uk](http://www.venhealth.co.uk)

 **Laboratoire Unither** 151 rue André Durouchez, Espace Industriel Nord,  
CS 28028, 80 084 Amiens Cedex 2 - France

Version of instructions  
for use : 02/2014



**VENTURE**  
HEALTHCARE<sup>LTD</sup>





### Indications:

Increases the mobilisation of secretions in the lower respiratory tract by osmotic effects. Prevents drying of bronchial mucous. This hypertonic saline solution 6% is for inhalation and can be used with pneumatic or ultrasonic nebulisers. Use your nebuliser according to the manufacturer's instructions provided with the system.

### Forms and presentations:

Pack of 60 single-dose vials of 4mL.

### Composition:

Sodium chloride: 6g, purified water qs 100mL. Preservative free.

### Directions for use:

The recommended dose of Hypertonic saline solution 6% for inhalation is one unit-dose vial 2 to 4 times a day. Break off a unit-dose vial and open it by twisting off the top. Pour the inhalation solution into the nebuliser. Use your nebuliser according to the manufacturer's instructions provided with the system.

### Side effects:

In rare cases, temporary irritations (e.g. coughing) or reversible constriction of the bronchia can occur. Pretreatment with bronchodilatory medication is recommended. If undesirable side effects do occur, discontinue treatment and seek advice from your doctor or pharmacist.

### Precaution for use:

Single use only. For inhalation only. Not for injection. Check that the unit-dose vial has not been damaged prior to use. Do not use a damaged unit-dose vial. Follow the nebuliser manufacturer's instructions for use. If a tendency to develop dyspnea or hypersensitivity is present, the inhalation solution must not be used without medical supervision. Do not mix with other medications. Keep out of the reach of children. Do not use after the expiry date. Opened single-dose ampoules must not be reused: since the product does not contain preservatives a re-use would lead to the risk of administration of a contaminated solution.

The hypertonic saline solution 6% for inhalation is a Class IIa medical device. Date of CE marking: 27/02/2014.

#### Explanation of symbols:

Symbol for sterility.



Symbol for single use only.



Please follow the instructions for use.



Use by date.



Manufacturer.